

June 1, 2016



DR 4.3 Features:

- **Wide center rib:** for stability and even wear
- **Open shoulder:** provides good evacuation in adverse conditions
- **Open shoulder:** provides optimal traction and grip for chains
- **Trusted, proprietary compound:** proven performance.
- **New, modern design:** balances price and functionality

DR 4.3

New Product

Ordering Information:

Tread Size	Item Number	Skid Quantity	Tread Depth
190mm	3001689	20	22/32
200mm	3001690	20	22/32
210mm	3001691	20	22/32
220mm	3001692	20	22/32
230mm	3001693	20	22/32
240mm	3001694	20	22/32
250mm	3001695	20	22/32

Curing Information:

Rim, Flange, Tube and Envelope	Inside Envelope	ARC™
3 hrs. 50 min.	3 hrs. 50 min.	3 hrs. 30 min.

Other Information:

OPTI-CURE™	Strip Wick	QuikCure™
Standard	Yes	Yes

For Radial Tires Only.

DR 4.3 is eligible for Smart Resource and Q-Fund.